

Created September 2016.

REF 07P8722

REF 07P8732

Instructions must be carefully followed. Reliability of assay results cannot be guaranteed if there are any deviations from these instructions.

### Name

Alinity i Anti-HBc II Reagent Kit (also referred to as Anti-HBc)

### INTENDED USE

The Alinity i Anti-HBc II assay is a chemiluminescent microparticle immunoassay (CMIA) used for the qualitative detection of antibody to hepatitis B core antigen (anti-HBc) in human serum and plasma, including specimens collected post-mortem (non-heartbeating) on the Alinity i analyzer.

The Alinity i Anti-HBc II assay is intended to be used as an aid in the diagnosis of hepatitis B infection and as a screening test to prevent transmission of hepatitis B virus (HBV) to recipients of blood, blood components, cells, tissue and organs.

### SUMMARY AND EXPLANATION OF THE TEST

The Alinity i Anti-HBc II assay utilizes microparticles coated with recombinant hepatitis B virus core antigen (rHBcAg) for the detection of anti-HBc. Anti-HBc determinations can be used as an indicator of current or past HBV infection. Anti-HBc is found in serum shortly after the appearance of hepatitis B surface antigen (HBsAg) in acute HBV infections. It will persist after the disappearance of HBsAg and before the appearance of detectable antibody to HBsAg (anti-HBs).<sup>1-7</sup> In the absence of information about any other HBV markers, it must be considered that an individual with detectable levels of anti-HBc may be actively infected with HBV or that the infection may have resolved, leaving the person immune.<sup>8</sup> Anti-HBc may be the only serological marker of HBV infection and potentially infectious blood.<sup>9-15</sup>

The presence of anti-HBc does not differentiate between acute or chronic hepatitis B infection.

### BIOLOGICAL PRINCIPLES OF THE PROCEDURE

This assay is a two-step immunoassay for the qualitative detection of anti-HBc in human serum and plasma using chemiluminescent microparticle immunoassay (CMIA) technology.

Sample, rHBcAg coated paramagnetic microparticles, specimen diluent, and assay diluent are combined and incubated. The anti-HBc present in the sample binds to the rHBcAg coated microparticles. The mixture is washed. Anti-human acridinium-labeled conjugate is added to create a reaction mixture and incubated. Following a wash cycle, Pre-Trigger and Trigger Solutions are added.

The resulting chemiluminescent reaction is measured as relative light units (RLUs). There is a direct relationship between the amount of anti-HBc in the sample and the RLUs detected by the system optics. The presence or absence of anti-HBc in the sample is determined by comparing the chemiluminescent RLU in the reaction to the cutoff RLU determined from an active calibration.

If the chemiluminescent signal in the reaction is greater than or equal to the cutoff signal, the specimen is considered reactive for anti-HBc.

For additional information on system and assay technology, refer to the Alinity ci-series Operations Manual, Section 3.

### REAGENTS

#### Kit Contents

Alinity i Anti-HBc II Reagent Kit 07P87

NOTE: Some kit sizes are not available in all countries. Please contact your local distributor.

NOTE: This product is composed of 4 components, which are packaged as a 2 cartridge reagent set. Both cartridges are required to perform the assay.

Volumes (mL) listed in the table below indicate the volume per cartridge set.


REF	07P8722	07P8732
Tests per cartridge set	100	600
Number of cartridge sets per kit	2	2
Tests per kit	200	1200
<b>MICROPARTICLES</b>	4.2 mL	16.8 mL
<b>CONJUGATE</b>	8.7 mL	17.8 mL
<b>ASSAY DILUENT</b>	5.9 mL	15.3 mL
<b>SPECIMEN DILUENT</b>	4.2 mL	15.1 mL
<b>MICROPARTICLES</b> Hepatitis B core ( <i>E. coli</i> , recombinant) antigen coated microparticles in TRIS buffer. Minimum concentration: 0.08% solids. Preservatives: ProClin 950 and sodium azide.		
<b>CONJUGATE</b> Murine acridinium-labeled anti-human conjugate in MES buffer with protein stabilizers. Minimum concentration: 0.04 µg/mL. Preservatives: sodium alkyl paraben and sodium azide.		
<b>ASSAY DILUENT</b> Murine protein stabilizers in MOPSO buffer. Preservatives: ProClin 950 and sodium azide.		
<b>SPECIMEN DILUENT</b> Reductant in MOPSO buffer.		

#### Warnings and Precautions



- **IVD**
- For *In Vitro* Diagnostic Use

#### Safety Precautions

**CAUTION:** This product requires the handling of human specimens. It is recommended that all human-sourced materials be considered potentially infectious and handled in accordance with the OSHA Standard on Bloodborne Pathogens. Biosafety Level 2 or other appropriate biosafety practices should be used for materials that contain or are suspected of containing infectious agents.<sup>16-19</sup>

The following warnings and precautions apply to: <b>MICROPARTICLES</b>	
	
<b>WARNING</b>	Contains methylisothiazolone and sodium azide.
H317	May cause an allergic skin reaction.
EUH032	Contact with acids liberates very toxic gas.
<b>Prevention</b>	
P261	Avoid breathing mist / vapors / spray.
P272	Contaminated work clothing should not be allowed out of the workplace.
P280	Wear protective gloves / protective clothing / eye protection.
<b>Response</b>	
P302+P352	IF ON SKIN: Wash with plenty of water.
P333+P313	If skin irritation or rash occurs: Get medical advice / attention.
P362+P364	Take off contaminated clothing and wash it before reuse.
<b>Disposal</b>	
P501	Dispose of contents / container in accordance with local regulations.

The following warnings and precautions apply to: <b>CONJUGATE</b>	
Contains sodium azide.	
EUH032	Contact with acids liberates very toxic gas.
P501	Dispose of contents / container in accordance with local regulations.

The following warnings and precautions apply to: <b>ASSAY DILUENT</b>	
 	
<b>DANGER</b>	Contains polyethylene glycol octylphenyl ether (Triton X-405), methylisothiazolone and sodium azide.
H318	Causes serious eye damage.
H317	May cause an allergic skin reaction.
H412	Harmful to aquatic life with long lasting effects.
EUH032	Contact with acids liberates very toxic gas.
<b>Prevention</b>	
P261	Avoid breathing mist / vapors / spray.
P272	Contaminated work clothing should not be allowed out of the workplace.
P280	Wear protective gloves / protective clothing / eye protection.
P273	Avoid release to the environment.
<b>Response</b>	
P302+P352	IF ON SKIN: Wash with plenty of water.
P305+P351+P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P310	Immediately call a POISON CENTER or doctor / physician.
P333+P313	If skin irritation or rash occurs: Get medical advice / attention.
P362+P364	Take off contaminated clothing and wash it before reuse.
<b>Disposal</b>	
P501	Dispose of contents / container in accordance with local regulations.

Safety Data Sheets are available at [www.abbottdiagnostics.com](http://www.abbottdiagnostics.com) or contact your local representative.

For a detailed discussion of safety precautions during system operation, refer to the Alinity ci-series Operations Manual, Section 8.

#### Reagent Handling

- Upon receipt, gently invert the unopened reagent kit by rotating it over and back for a full 180 degrees, 5 times with green label stripe facing up and then 5 times with green label stripe facing down. This ensures that liquid covers all sides of the bottles within the cartridges. During reagent shipment, microparticles can settle on the reagent septum.
  - Place a check in the square on the reagent kit to indicate to others that the inversions have been completed.
- After mixing, place reagent cartridges in an upright position for 1 hour before use to allow bubbles that may have formed to dissipate.
- If a reagent cartridge is dropped, place in an upright position for 1 hour before use to allow bubbles that may have formed to dissipate.
- Reagents are susceptible to the formation of foam and bubbles. Bubbles may interfere with the detection of the reagent level in the cartridge and cause insufficient reagent aspiration that may adversely affect results.

For a detailed discussion of reagent handling precautions during system operation, refer to the Alinity ci-series Operations Manual, Section 7.

#### Reagent Storage

	Storage Temperature	Maximum Storage Time	Additional Storage Instructions
<b>Unopened</b>	2 to 8°C	Until expiration date	Store in upright position. If cartridge does not remain upright, gently invert the cartridge 10 times and place in an upright position for 1 hour before use.
<b>Onboard</b>	System Temperature	30 days	
<b>Opened</b>	2 to 8°C	Until expiration date	Store in upright position. If cartridge does not remain upright during storage, discard the cartridge. Do not reuse original reagent caps or replacement caps due to the risk of contamination and the potential to compromise reagent performance.

Reagents may be stored on or off the system. If removed from the system, store reagents with new replacement caps in an upright position at 2 to 8°C. For reagents stored off the system, it is recommended that they be stored in their original trays or boxes to ensure they remain upright.

For information on unloading reagents, refer to the Alinity ci-series Operations Manual, Section 5.

#### Indications of Reagent Deterioration

Deterioration of the reagents may be indicated when a calibration error occurs or a control value is out of the specified range. Associated test results are invalid, and samples must be retested. Assay recalibration may be necessary.

For troubleshooting information, refer to the Alinity ci-series Operations Manual, Section 10.

## INSTRUMENT PROCEDURE

The Alinity i Anti-HBc II assay file must be installed on the Alinity i analyzer prior to performing the assay.

For detailed information on assay file installation and viewing and editing assay parameters, refer to the Alinity ci-series Operations Manual, Section 2.

For information on printing assay parameters, refer to the Alinity ci-series Operations Manual, Section 5.

For a detailed description of system procedures, refer to the Alinity ci-series Operations Manual.

## SPECIMEN COLLECTION AND PREPARATION FOR ANALYSIS

### Specimen Types

The specimen types listed below were verified for use with this assay on the ARCHITECT i System.

Other specimen types and collection tube types have not been verified with this assay.

Specimen Types	Collection Tubes
Serum	Serum
	Serum separator
Plasma	Sodium heparin
	Lithium heparin (PST)
	Potassium-EDTA
	Sodium citrate
	Potassium oxalate
	CPD
	CPDA-1
	ACD

- ACD tubes may show a positive bias up to 20% relative to serum.
- Liquid anticoagulants may have a dilution effect resulting in lower concentrations for individual patient specimens.
- Performance has been established for the use of cadaveric blood specimens (specimens collected post-mortem, non-heart-beating) that have been collected up to 17.5 hours after death. Performance was established using 50 spiked and 50 non-spiked cadaveric blood specimens.<sup>20</sup>
- Testing of cadaveric blood specimens from patients with plasma dilution due to transfusions of > 2000 mL of blood or colloids within 48 hours, or > 2000 mL of crystalloids within 1 hour (or any combination thereof) prior to collection of the specimens have not been validated.
- For cadaveric blood specimens follow general standards and/or regulations for collection, storage and handling.
- The instrument does not provide the capability to verify specimen types. It is the responsibility of the operator to verify that the correct specimen types are used in the assay.

### Specimen Conditions

- Do not use:
  - heat-inactivated specimens
  - pooled specimens
  - grossly hemolyzed specimens
  - specimens with obvious microbial contamination
  - body fluids other than human serum and plasma
- For accurate results, serum and plasma specimens should be free of fibrin, red blood cells, and other particulate matter. Serum specimens from patients receiving anticoagulant or thrombolytic therapy may contain fibrin due to incomplete clot formation.
- To prevent cross contamination, use of disposable pipettes or pipette tips is recommended.

### Preparation for Analysis

- Follow the tube manufacturer's processing instructions for collection tubes. Gravity separation is not sufficient for specimen preparation.
- Specimens should be free of bubbles. Remove bubbles with an applicator stick before analysis. Use a new applicator stick for each specimen to prevent cross-contamination.

To ensure consistency in results, recentrifuge specimens prior to testing if

- they contain fibrin, red blood cells, or other particulate matter
- they require repeat testing.

NOTE: If fibrin, red blood cells, or other particulate matter are observed, mix by low speed vortex or by inverting 10 times prior to recentrifugation.

Prepare frozen specimens as follows:

- Frozen specimens must be completely thawed before mixing.
- Mix thawed specimens thoroughly by low speed vortexing or by inverting 10 times.
- Visually inspect the specimens. If layering or stratification is observed, mix until specimens are visibly homogeneous.
- If specimens are not mixed thoroughly, inconsistent results may be obtained.
- Recentrifuge specimens.

Prepare cadaveric blood specimens as follows:

- After initial centrifugation, recentrifuge specimens as described below.
- If specimens are not processed directly after initial centrifugation, it is recommended to remove the supernatant from the clot, red blood cells or separator gel until further processing.

### Re-centrifugation of Specimens

- Transfer mixed specimens to a centrifuge tube and centrifuge at a minimum of 100 000 g-minutes.
- Examples of acceptable time and force ranges that meet this criterion are listed in the table below.

Centrifugation time using alternate RCF values can be calculated using the following formula:

$$\text{Minimum Centrifugation time (minutes)} = \frac{100\,000 \text{ g-minutes}}{\text{RCF}}$$

Recentrifugation Time (Minutes)	RCF (x g)	g-Minutes
10	10 000	100 000
20	5000	100 000
40	2500	100 000

$$\text{RCF} = 1.12 \times r_{\text{max}} (\text{rpm}/1000)^2$$

RCF - The relative centrifugal force generated during centrifugation.

rpm - The revolutions per minute of the rotor on which the specimens are being spun (usually the digital readout on the centrifuge will indicate the rpm).

Centrifugation Time - The time should be measured from the time the rotor reaches the required RCF or rpm to the time it begins decelerating.

$r_{\text{max}}$  - Radius of the rotor in millimeters. **NOTE:** If custom tube adapters (i.e., adapters not defined by the centrifuge manufacturer) are used, then the radius (r<sub>max</sub>) should be manually measured in millimeters and the RCF calculated.

g-minutes - The unit of measure for the product of RCF (x g) and centrifugation time (minutes).

- Transfer clarified specimen to a sample cup or secondary tube for testing. For centrifuged specimens with a lipid layer, transfer only the clarified specimen and not the lipemic material.

## Specimen Storage

Specimen Type	Temperature	Maximum Storage Time	Special Instructions
Serum/ Plasma	Room temperature (15 to 30°C)	3 days	Specimens may be stored on or off the clot, red blood cells, or separator gel.
	2 to 8°C	14 days	Specimens may be stored on or off the clot, red blood cells, or separator gel.
Cadaveric	Room temperature (15 to 30°C)	3 days	If specimens are not processed directly after initial centrifugation, it is recommended to remove the supernatant from the clot, red blood cells or separator gel until further processing.
	2 to 8°C	7 days	If specimens are not processed directly after initial centrifugation, it is recommended to remove the supernatant from the clot, red blood cells or separator gel until further processing.

Remove serum or plasma from the clot, red blood cells, or separator gel if stored longer than the maximum 15-30°C or 2-8°C storage time and store frozen at -20°C or colder.

For serum and plasma no qualitative performance differences were observed between experimental controls and nonreactive or spiked reactive specimens subjected to 6 freeze/thaw cycles.

No qualitative differences were observed for cadaveric blood specimens (nonreactive or spiked reactive) when subjected to up to 3 freeze/thaw cycles.

Avoid multiple freeze/thaw cycles for all specimen types.

## Specimen Shipping

Package and label specimens in compliance with applicable state, federal, and international regulations covering the transport of clinical specimens and infectious substances.

## PROCEDURE

### Materials Provided

07P87 Alinity i Anti-HBc II Reagent Kit

### Materials Required but not Provided

- Alinity i Anti-HBc II assay file
- 07P8701 Alinity i Anti-HBc II Calibrator
- 07P8710 Alinity i Anti-HBc II Controls or other control material
- 06P1160 Alinity Trigger Solution
- 06P1265 Alinity Pre-Trigger Solution
- 06P1368 Alinity i-series Concentrated Wash Buffer

For information on materials required for operation of the instrument, refer to the Alinity ci-series Operations Manual, Section 1.

For information on materials required for maintenance procedures, refer to the Alinity ci-series Operations Manual, Section 9.

## Assay Procedure

For a detailed description of how to run an assay, refer to the Alinity ci-series Operations Manual, Section 5.

- If using primary or aliquot tubes, refer to the Alinity ci-series Operations Manual, Section 4 to ensure sufficient specimen is present.
- To minimize the effects of evaporation, verify adequate sample cup volume is present prior to running the test.
- Maximum number of replicates sampled from the same sample cup: 10
  - Priority:
    - Sample volume for first test: 75 µL
    - Sample volume for each additional test from same sample cup: 25 µL
  - ≤ 3 hours on the reagent and sample manager:
    - Sample volume for first test: 150 µL
    - Sample volume for each additional test from same sample cup: 25 µL
  - > 3 hours on the reagent and sample manager:
    - Replace with a fresh aliquot of sample.
- Refer to the Alinity i Anti-HBc II calibrator package insert and Alinity i Anti-HBc II control package insert for preparation and usage.
- For general operating procedures, refer to the Alinity ci-series Operations Manual, Section 5.
- For optimal performance, it is important to perform routine maintenance as described in the Alinity ci-series Operations Manual, Section 9. Perform maintenance more frequently when required by laboratory procedures.

## Sample Dilution Procedures

Samples cannot be diluted for the Alinity i Anti-HBc II assay.

## Calibration

For instructions on performing a calibration, refer to the Alinity ci-series Operations Manual, Section 5.

Each assay control must be tested to evaluate the assay calibration.

Once a calibration is accepted and stored, all subsequent samples may be tested without further calibration unless:

- A reagent kit with a new lot number is used.
- Daily quality control results are outside of statistically-based quality control limits used to monitor and control system performance, as described in the Quality Control Procedures section of this package insert.
  - If statistically-based quality control limits are not available then the calibration should not exceed a 30-day limit for recalibration frequency.

This assay may require recalibration after maintenance to critical parts or subsystems or after service procedures have been performed.

## Quality Control Procedures

The recommended control requirement for the Alinity i Anti-HBc II assay is that a single sample of each control level be tested once every 24 hours each day of use.

Additional controls may be tested in accordance with local, state, and/or federal regulations or accreditation requirements and your laboratory's quality control policy.

To establish statistically-based control limits, each laboratory should establish its own concentration target and ranges for new control lots at each clinically relevant control level. This can be accomplished by assaying a minimum of 20 replicates over several (3-5) days and using the reported results to establish the expected average (target) and variability about this average (range) for the laboratory. Sources of variation that should be included in this study in order to be representative of future system performance include:

- Multiple stored calibrations
- Multiple reagent lots
- Multiple calibrator lots
- Multiple processing modules (if applicable)
- Data points collected at different times of the day

Refer to published guidelines for information or general control recommendation, for example Clinical and Laboratory Standards Institute (CLSI) Document C24-A3 or other published guidelines, for general quality control recommendations.<sup>21</sup>

- If more frequent control monitoring is required, follow the established quality control procedures for your laboratory.
- If quality control results do not meet the acceptance criteria defined by your laboratory, sample results may be suspect. Follow the established quality control procedures for your laboratory. Recalibration may be necessary. For troubleshooting information, refer to the Alinity ci-series Operations Manual, Section 10.
- Review quality control results and acceptance criteria following a change of reagent or calibrator lot.

#### Quality Control Guidance

Refer to “Basic QC Practices” by James O Westgard, Ph.D. for guidance on laboratory quality control practices.<sup>22</sup>

#### Verification of Assay Claims

For protocols to verify package insert claims, refer to Verification of Assay Claims in the Alinity ci-series Operations Manual.

## RESULTS

#### Calculation

The Alinity i analyzer calculates results for the Alinity i Anti-HBc II assay using the ratio of the sample RLU to the cutoff RLU (S/CO) for each specimen and control.

Cutoff RLU = Calibrator 1 Mean RLU x 1.0.

The cutoff RLU is stored for each reagent lot calibration.

S/CO = Sample RLU/Cutoff RLU

#### Interpretation of Results

The cutoff is 1.00 S/CO.

Initial Results		
S/CO	Instrument Interpretation	Retest Procedure
< 1.00	Nonreactive	No retest required.
≥ 1.00	Reactive	Retest in duplicate.
Final Interpretation		
Initial Interpretation	Results with Retest	Final Interpretation
Nonreactive	No retest required.	Nonreactive
Reactive	If two of the three results are < 1.00 S/CO	Nonreactive
	If two of the three results are ≥ 1.00 S/CO	Reactive

For details on configuring the Alinity i analyzer to use grayzone interpretations, refer to the Alinity ci-series System Operations Manual.

The grayzone interpretation from the Alinity i interpretations screen is not used by the Alinity i analyzer unless a grayzone is configured.

#### Flags

Some results may contain information in the Flags field. For a description of the flags that may appear in this field, refer to the Alinity ci-series Operations Manual, Section 5.

## LIMITATIONS OF THE PROCEDURE

- If the anti-HBc results are inconsistent with clinical evidence, additional testing is suggested to confirm the result.
- For diagnostic purposes, results should be used in conjunction with other data; e.g., symptoms, results of other tests, clinical impressions, etc.
- Heterophilic antibodies in human serum can react with reagent immunoglobulins, interfering with in vitro immunoassays. Patients routinely exposed to animals or to animal serum products can be prone to this interference and anomalous values may be observed. Additional information may be required for diagnosis.<sup>23</sup>
- Specimens from patients who have received preparations of mouse monoclonal antibodies for diagnosis or therapy may contain human anti-mouse antibodies (HAMA). Specimens containing HAMA may produce anomalous values when tested with assay kits that employ mouse monoclonal antibodies.<sup>24, 25</sup>

## SPECIFIC PERFORMANCE CHARACTERISTICS

Representative performance data are provided in this section. Results obtained in individual laboratories may vary.

The Alinity i analyzer and the ARCHITECT i System utilize the same reagents and sample/reagent ratios.

Unless otherwise specified, all studies were performed on the Alinity i analyzer.

#### Precision

##### Within-Laboratory Precision

A study was performed based on guidance from CLSI EP05-A2.<sup>26</sup>

Testing was conducted using 3 lots of the Alinity i Anti-HBc II Reagent Kit, 3 lots of the Alinity i Anti-HBc II Calibrator, and 3 lots of the Alinity i Anti-HBc II Controls and 1 instrument. Two controls and 2 human plasma panels were assayed in a minimum of 2 replicates at 2 separate times per day on 20 different days.

Sample	n	Mean (S/CO)	Within-Run (Repeatability)		Within-Laboratory (Total) <sup>a</sup>	
			SD	%CV	SD (Range <sup>b</sup> )	%CV (Range <sup>b</sup> )
Negative Control	359	0.12	0.009	7.2	0.010 (0.008-0.012)	8.3 (6.2-9.5)
Positive Control	357	2.78	0.062	2.2	0.095 (0.092-0.098)	3.4 (3.3-3.5)
Panel 1	350	0.83	0.020	2.4	0.029 (0.027-0.034)	3.6 (3.3-3.9)
Panel 2	358	1.23	0.032	2.6	0.045 (0.042-0.049)	3.6 (3.4-3.9)

<sup>a</sup> Includes within-run, between-run, and between-day variability.

<sup>b</sup> Maximum and minimum SD or %CV for each reagent lot and instrument combination.

#### Specificity

A total of 5391 specimens, including specimens from blood donors and hospitalized patients were tested using the Alinity i Anti-HBc II assay and a commercially-available anti-HBc assay. Repeatedly reactive samples were further tested by supplemental testing.



Category	Alinity i Anti-HBc II					Commercially-Available Anti-HBc Assay	
	n	Number		Specificity <sup>a</sup> (95% CI)	Specificity <sup>a</sup> (95% CI)	n	Specificity <sup>a</sup> (95% CI)
		IR (% of Total)	RR (% of Total)				
Blood Donors - Serum	2625	6 (0.23)	4 (0.15)	0 (2621/2625) (99.61 - 99.96)	99.85% (2625/2627) (99.73 - 99.99)	2627	99.92% (2625/2627) (99.73 - 99.99)
Blood Donors - Plasma	2548	7 (0.27)	7 (0.27)	4 (2541/2544) (99.66 - 99.98)	99.88% (2540/2544) (99.60 - 99.96)	2548	99.84% (2540/2544) (99.60 - 99.96)
Total Donors	5173	13 (0.25)	11 (0.21)	4 (5162/5169) (99.72 - 99.95)	99.86% (5165/5171) (99.75 - 99.96)	5175	99.88% (5165/5171) (99.75 - 99.96)
Hospitalized Patients	218	7 (3.21)	7 (3.21)	7 (211/211) (98.27 - 100.00)	100.00% (211/211) (98.27 - 100.00)	218	100.00% (211/211) (98.27 - 100.00)

IR = Initially Reactive, RR = Repeatedly Reactive, CI = Confidence Interval

<sup>a</sup> Repeatedly reactive specimens determined to be positive by supplemental testing were excluded from these calculations.

### Sensitivity

A total of 408 anti-HBc positive specimens from patients with acute, chronic and past/resolved HBV infection were tested using the Alinity i Anti-HBc II assay and a commercially-available anti-HBc assay, resulting in a sensitivity of 100.00% (408/408) which is within the 95% confidence interval of 99.10% to 100% with a 100% expected sensitivity.

Specimen Category	n	Number	RR	Alinity i Anti-HBc II Sensitivity	Commercially-Available Anti-HBc Assay Sensitivity
Anti-HBc Positive	408	408		100.00%	100.00%

RR = Repeatedly Reactive

### Analytical Sensitivity

Analytical sensitivity was evaluated using dilutions of the WHO 1st International Standard for anti-Hepatitis B core antigen (anti-HBc), plasma, human NIBSC code: 95/522. The dilutions ranged from 1.00 to 0.05 IU/mL. The dilutions were tested across 3 reagent lots on 1 instrument. The analytical sensitivity results on the Alinity i Anti-HBc II assay ranged from 0.54 to 0.56 IU/mL.

Note: The 1st WHO International Standard for anti-Hepatitis B core antigen (anti-HBc), NIBSC 95/522, was assessed against the Paul-Ehrlich-Institute (PEI) anti-HBc standard (PEI 82). PEI units (PEI U) were found to be equivalent to the International Units (IU).<sup>27</sup>

### Other Disease States

These studies were performed on the ARCHITECT i System.

Additional studies were performed to evaluate other potential interfering disease states on the ARCHITECT Anti-HBc II assay. A total of 104 specimens were tested from the following categories: antinuclear antibodies (ANA), Epstein-Barr virus (anti-EBV positive), hepatitis A virus (anti-HAV IgM positive), hepatitis C virus (anti-HCV positive), human immunodeficiency virus (anti-HIV-1 positive), human anti-mouse antibodies (HAMA) positive, influenza vaccine recipients, non-viral liver disease, rheumatoid factor positive, syphilis, systemic lupus erythematosus (SLE), toxoplasmosis IgG positive, varicella zoster (anti-VZV positive), anti-*E. coli* positive and yeast infection. With these specimens, ARCHITECT Anti-HBc II showed the same qualitative results as the comparator method.

### Interference

This study was performed on the ARCHITECT i System.

#### Potentially Interfering Endogenous Substances

No qualitative performance differences were observed between experimental controls and nonreactive or reactive specimens tested with the compounds listed in the table below.

Potentially Interfering Substance	Interferent Level
Bilirubin	≤ 20 mg/dL
Triglycerides	≤ 3000 mg/dL
Protein	≤ 12 g/dL
Hemoglobin	≤ 500 mg/dL

### BIBLIOGRAPHY






1. Hoofnagle JH, Gerety RJ, Barker LF. Antibody to hepatitis-B-virus core in man. *Lancet* 1973;ii:869-873.
2. Szmunes W, Hoofnagle JH, Stevens CE, et al. Antibody against the hepatitis type B core antigen. A new tool for epidemiologic studies. *Am J Epidemiol* 1976;104(3):256-262.
3. Hoofnagle JH, Seeff LB, Buskell-Bales Z, et al. Serologic responses in HB. In: Vyas GN, Cohen SN, Schmid R, editors. *Viral Hepatitis*. Philadelphia, PA: Franklin Institute Press; 1978:219-242.
4. Krugman S, Overby LR, Mushahwar IK, et al. Viral hepatitis, type B studies on natural history and prevention re-examined. *N Engl J Med* 1979;300:101-106.
5. Zito DR, Gurdak RG, Tucker FL, et al. Hepatitis B virus serology: Loss of antibody to surface antigen. *Am J Clin Pathol* 1987;88(2):229-231.
6. Gitlin N. Hepatitis B; Diagnosis, Prevention, and Treatment. *Clin Chem* 1997;43(8B):1500-1506.
7. Koff RS. Viral Hepatitis. In: Schiff L, Schiff ER, editors. *Diseases of the Liver*. 7th ed. Philadelphia: JB Lippincott, 1993:492-577.
8. Dodd RY, Popovsky MA, Members of the Scientific Section Coordinating Committee. Antibodies to hepatitis B core antigen and the infectivity of the blood supply. *Transfusion* 1991;31(5):443-449.
9. Seeff LB, Beebe GW, Hoofnagle JH, et al. A serologic follow-up of the 1942 epidemic of post-vaccination hepatitis in the United States Army. *N Engl J Med* 1987;316(16):965-970.
10. Katchaki JN, Siem TH, Brouwer R, et al. Detection and significance of anti-HBc in the blood bank; preliminary results of a controlled prospective study. *J Virol Methods* 1980;2:119-125.
11. Lander JJ, Gitnick GL, Gelb LH, et al. Anticore antibody screening of transfused blood. *Vox Sang* 1978;34:77-80.
12. Hoofnagle JH, Seeff LB, Buskell-Bales Z, et al. Type B hepatitis after transfusion with blood containing antibody to hepatitis B core antigen. *N Engl J Med* 1978;298(25):1379-1383.
13. Koziol DE, Holland PV, Alling DW, et al. Antibody to hepatitis B core antigen as a paradoxical marker for non-A, non-B hepatitis agents in donated blood. *Ann Intern Med* 1986;104:488-495.
14. AuBuchon JP, Sandler SG, Fang CT, et al. American Red Cross experience with routine testing for hepatitis B core antibody. *Transfusion* 1989;29:230-232.
15. Stevens CE, Aach RD, Hollinger FB, et al. Hepatitis B virus antibody in blood donors and the occurrence of non-A, non-B hepatitis in transfusion recipients: An analysis of the transfusion-transmitted viruses study. *Ann Intern Med* 1984;101(6):733-738.
16. US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
17. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
18. World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
19. Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.
20. U.S. Department of Health and Human Services, Food and Drug Administration, Center for Biologics Evaluation and Research. Guidance for Industry Recommendations for Obtaining a Labeling Claim for Communicable Disease Donor Screening Tests Using Cadaveric Blood Specimens from Donors of Human Cells, Tissues, and Cellular and Tissue-Based Products (HCT/PS), November 2004. <http://www.fda.gov/BiologicsBloodVaccines/GuidanceComplianceRegulatoryInformation/Guidances/Tissue/ucm073972.htm> Accessed June 1, 2016

21. Clinical and Laboratory Standards Institute (CLSI). *Statistical Quality Control for Quantitative Measurement Procedures: Principles and Definitions; Approved Guideline—Third Edition*. CLSI Document C24-A3. Wayne, PA: CLSI; 2006.
22. Westgard JO. *Basic QC Practices*. 3rd ed. Madison, WI: Westgard Quality Corporation; 2010.
23. Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem* 1988;34(1):27-33.
24. Primus FJ, Kelley EA, Hansen HJ, et al. "Sandwich"-type immunoassay of carcinoembryonic antigen in patients receiving murine monoclonal antibodies for diagnosis and therapy. *Clin Chem* 1988;34(2):261-264.
25. Schroff RW, Foon KA, Beatty SM, et al. Human anti-murine immunoglobulin responses in patients receiving monoclonal antibody therapy. *Cancer Res* 1985;45(2):879-885.
26. Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition*. CLSI Document EP05-A2. Wayne, PA: CLSI; 2004.
27. Scheiblaue H, Ferguson M, Volkers P, Nick S. Report of the WHO collaborative study to establish the First International Standard for detection of antibodies to Hepatitis B core antigen (anti-HBc), human plasma. WHO Report, WHO/BS/08.2098, October 2008.

Note for number formatting:

- A space is used as thousands separator (example: 10 000 specimens).
- A period is used to separate the integer part from the fractional part of a number written in decimal form (example: 3.12%).

## ■ Key to Symbols

	Consult instructions for use
	Manufacturer
	Sufficient for
	Temperature limitation
	Use by/Expiration date
<b>ASSAY DILUENT</b>	Assay Diluent
<b>CONJUGATE</b>	Conjugate
<b>CONTAINS: AZIDE</b>	Contains Sodium Azide. Contact with acids liberates very toxic gas.
<b>INVERSIONS PERFORMED</b>	Inversions Performed
<b>IVD</b>	<i>In Vitro</i> Diagnostic Medical Device
<b>LOT</b>	Lot Number
<b>MICROPARTICLES</b>	Microparticles
<b>PRODUCT OF GERMANY</b>	Product of Germany
<b>REF</b>	List Number
<b>SN</b>	Serial number
<b>SPECIMEN DILUENT</b>	Specimen Diluent

Alinity and ARCHITECT are trademarks of Abbott Laboratories in various jurisdictions. All other trademarks are property of their respective owners.



Abbott GmbH & Co. KG  
Max-Planck-Ring 2  
65205 Wiesbaden  
Germany  
+49-6122-580



0088

**Customer Service: Contact your local representative  
or find country-specific contact information on  
[www.abbottdiagnostics.com](http://www.abbottdiagnostics.com)**

Created September 2016.

©2016 Abbott Laboratories